



Release Notes

Medtech32

Web Form Validation Update

(January 2012)



IMPORTANT NOTE

Medtech recommends that all Medtech upgrades and database back-up and restore processes are performed by a **Medtech Certified Technician** or your Practice Technician.

For a list of Medtech Certified technicians please visit the Medtech Global Website:

www.medtechglobal.com

These Release Notes contain important information for all Medtech Users. Please ensure that they are circulated amongst all your staff. We suggest that these notes are filed safely for future reference.

Issues Resolved in this Update

1. Web Form Validation

This update resolves the “advanced form validation” error encountered by some practices when attempting to run the following Web Forms.

- ADHB Referral
- General Referral
- MSD Medical Certificate
- Hutt Valley DHB Referral
- HealthLink Online
- Konnect SureMed PMAR

Web Form Validation Update Installation

1. This update will run on any version of Medtech32.
2. Visit our website at www.medtechglobal.com.
3. Select “New Zealand” in the region drop down (top right hand side of the window).
4. Select menu, **Medtech Online** ▶ **Medtech32** ▶ **Downloads**. The New Zealand Medtech32 Updates page will be displayed.
5. Here you will find the **Web Form Validation Update**.
6. In the **File Download Security Warning** dialogue box, select **Run**.
7. In the **Open File Security Warning Dialogue** box, select **Run** option.
8. The installation screen will appear.
The logo for Medtech32, featuring the word "MEDtech" in a blue sans-serif font and "32" in a larger, bold blue font, all enclosed in a thin black rectangular border.
9. The installation process will begin.
10. Once the update is completed, the installation screen will close.

For further information on these new features, or any other queries regarding this update, please contact Medtech Support on 0800 2 MEDTECH (633 832) - Option 1 or send us a fax, toll free on 0800 MEDTECH (633 832) or email: support@medtechglobal.com.