



RELEASE NOTES

Medtech32

MeNZB/ACC Subsidy Increase & Clinical Event Export Utility Update

(March 2011)



IMPORTANT NOTE

This release **does not** require a Technician, however, if you choose to use a Technician, Medtech recommends that all Medtech upgrades and database back-up and restore processes are performed by a **Medtech Certified Technician** or your Practice Technician. For a list of Medtech Certified technicians please visit the Medtech Global Website: www.medtechglobal.com.

These Release Notes contain important information for all Medtech Users. Please ensure that they are circulated amongst all your staff. We suggest that these notes are filed safely for future reference.

Summary of Changes

1) ACC Rural Healthcare Subsidy Increase

There have been increases to the ACC Subsidy prices from 1st May 2011 for Rural Healthcare Providers. This update will provide users with the new subsidy prices with an effective date of 1st May 2011.

2) MeNZB Subsidy Removal

The Ministry of Health have advised that the MeNZB vaccine will no longer be available from 1 April 2011. Therefore no claim for the administration of the MeNZB vaccine from 1 April 2011 will be valid.

This update will set the subsidy for the administration of the MeNZB vaccine to \$0.00 with an effective date of 1st April 2011.

3) Clinical Event Export Utility

Currently there is a maximum 10-year rule applied to information extracted by the Medtech Clinical Performance Indicator Report this rule in some scenarios. This 10 year rule has now been removed.

Prerequisites

Please review the following prerequisites and ensure they are met prior to installation:

IMPORTANT: This documentation contains valuable information that, if not read, could seriously affect the upgrade progress and/or possible down time of your network.

- 1) This update will run on any version of Medtech32 greater than Version 18.0.
- 2) Ensure the person(s) who will be performing the upgrade have read through all of the release documentation and installation instructions.
- 3) Ensure you are logged onto Windows with Administrative rights when performing **ANY** installation, update, or maintenance tasks.
- 4) This update **MUST** be run on the Medtech32 Server. Users can remain logged into Medtech whilst this update is being run, however you must log out of the Medtech32 application and log in again to see the changes.

Medtech32 Update Installation

- 1) Ensure that you are running Version 18.0 or above of Medtech32. This can be checked by logging into Medtech32 and selecting **Help ► About Medtech32**.
- 2) Visit our website at www.medtechglobal.com.
- 3) Select 'New Zealand' in the region drop down (top right hand side of the window).
- 4) Select menu, **Medtech Online ► Medtech32 ► Downloads**. The New Zealand Medtech32 Updates page will be displayed.
- 5) Here you will find the **MeNZB/ACC Subsidy Increase & CPI Export Utility Update (March 2011)**
- 6) In the **File Download Security Warning** dialogue box, select **Run**.
- 7) In the **Open File Security Warning Dialogue** box, select **Run** option.
- 8) The installation process will begin.
- 9) Once the update is completed, the following message will be displayed.



Should you have any questions with regard to the information contained in these Release Notes please don't hesitate to contact the MedTech Helpdesk on **0800 2 MEDTECH (633 832) (Option 1 and follow prompts)**.