



Medtech32 Australia
ICPC-2 PLUS Update
(April 2012)

Release Notes



These Release Notes contain important information for all Medtech32 users.
Please ensure the Release Notes are circulated amongst all your staff.
We suggest these should be filed safely for future reference.

Table of Contents

For further information on this release, or any other queries regarding the ICPC-2 PLUS Update (April 2012), please contact the Medtech Helpdesk on 1800 148 165 → Option 1, or email ausupport@medtechglobal.com.

Table of Contents	2
Pre-requisites	3
Introduction	4
Installation	5
New Classification Codes	7

Pre-requisites

Please review the following pre-requisites and ensure they are met prior to running the Update:

- This Update can ONLY be run by **CURRENT** ICPC-2 PLUS Subscribers.
NOTE: For subscriptions enquiries, please contact the Family Medicine Research Centre at the University of Sydney on 02 9354 0616.

- Ensure the person(s) who will be performing the upgrade have **READ THROUGH** these notes.

IMPORTANT: This document contains valuable information that, if not read, could seriously affect the upgrade progress and/or possible down time of your network.

- Ensure you are ALWAYS logged onto Windows with **ADMINISTRATIVE RIGHTS** when performing ANY installation, update, or maintenance tasks.
- **"User Account Control" (UAC)** MUST be **DISABLED** as a policy across ALL Server and Client computers that are running on **Windows Vista or Windows 7 or Windows 2008**.
- To avoid unnecessary problems during installation, upgrade, or maintenance, preferably you should ALWAYS log onto Windows in **"Console Mode"** – i.e. not through Remote Desktop Connection, Terminal Services, or Citrix.
- Ensure you have a **COMPLETE** backup of ALL databases located in the MT32\Data directory.
- Ensure ALL users (including remote users) have **LOGGED OUT** of Medtech32 and ALL services (e.g. ManageMyHealth SMS Communicator), scheduled utilities (e.g. Message Transfer Utility, Auto Scan from Folder), backup or maintenance tasks (e.g. Interbase Backup Scripts) that require access to the databases have been **STOPPED**.
- Ensure the **LAST** Database you have logged in **LOCALLY** on the Medtech32 Server is the one you would like this update to run.
NOTE: If unsure, simply login to the CORRECT Database from the Medtech32 Server, then log out and close Medtech32 before continuing.

Introduction

This release notes provide an overview on how to apply the April 2012 ICPC-2 PLUS Update for Medtech32 ICPC-2 PLUS Subscribers.

IMPORTANT NOTE

WARNING: It is HIGHLY recommended to employ ONLY qualified system engineers when performing ANY installation and upgrade. The consequences of ruining a database during upgrade could possibly lead to data corruptions, and as a result – data loss and systems downtime.

If in doubt, please consult with your IT technician/service provider, or contact one of the Medtech Channel Partners listed on our web site:

<http://www.medtechglobal.com/aus/medtech-online-au/support-3.html>

IMPORTANT NOTE FOR BRIECASING

Changes from this Update will NOT be applied to the Briefcasing laptops UNTIL a **COMPLETE CHECK-OUT** has been run.

IMPORTANT: ALL Briefcasing laptops with OUTSTANDING Briefcased data MUST be checked-in PRIOR TO running a complete check-out.

Alternatively, if it is not feasible to run a complete check-out soon, you MUST install this Update separately on EACH Briefcasing laptop by following the same instructions in the "Installation" section below.

Installation

The ICPC-2 PLUS Update (April 2012) must be run on the Medtech32 Server machine. The following procedures ONLY need to be run ONCE for EACH practice (or once per database if your practice has multiple databases).

NOTE: If you are uncertain which computer is the Medtech32 Server, please contact your IT technician or service provider who has performed the Medtech32 installation and/or upgrade.


1. Please visit our website at www.medtechglobal.com.
2. Select **Region: Australia** from the Region dropdown list on the top right corner of the screen.




The Australia Home Page will be displayed.



3. Select from the Top Menu, **medtech online** ► **medtech32** ► **downloads** ► **third party database updates**. The Medtech32's Third Party Database Updates page will be displayed.



4. Here you will find the **ICPC2+ Update (April 2012)** and the **Release Notes**.
5. Click on the Link **ICPC2+ Update (April 2012)**.
6. If the File Download Security Warning dialogue box appears, select the **Run** option.
7. If the Open File Security Warning dialogue box appears, select the **Run** option.
8. The Update will then begin to run.



9. Once the Update is completed, the Update Complete screen will be displayed. Click on the Button **OK** to exit the Update.



New Classification Codes

The following is the list of new Classification Codes added by this ICPC-2 PLUS Update:

B79005	Accessory spleen
A85039	Adverse effect;anaesthetic
L45006	Advice/education;footwear
A45038	Advice/education;footwear
P45014	Advice/education;sleep
X99036	Anovulation
L99101	Arthritis;seronegative
U13009	Atonic;bladder
T02008	Binge eating
X99037	Blocked;fallopian tube
X21017	Calcification;breast
A13009	Concern (about);breast implant
A33059	Culture;gonococcal
L98018	Deformity;ankle;acquired
L82033	Deformity;ankle;congenital
K92031	Disease;small vessel
P73016	Disorder;affective
A59013	Electroacupuncture
D31007	Exam;rectal;digestive
D97013	Fibrosis;hepatic
W84030	Haematoma;subchorionic
L83030	Headache;cervicogenic
T93019	Hyperlipidaemia;familial
T99091	Hypoaldosteronism
D94011	Ileitis
A99025	Late effect;premature birth
D20038	Lesion;mouth
Z12029	Marriage;disruption
N89014	Migraine;cervicogenic
L57012	Mobilisation;joint
L02015	Pain;back;mechanical
L01011	Pain;neck;mechanical
K99053	Phlebectasis
S43003	Photography;skin
B99025	Polycythaemia
A94029	Problem;newborn
A94030	Problem;perinatal
D19033	Problem;wisdom teeth
R21019	Pus (on);tonsils
D67008	Referral;colorectal surgeon
S72003S	cabies;crusted
A85040	Side-effect;anaesthetic
S04018S	welling;groin
N94026	Tarsal tunnel syndrome

L92018	Tendonitis;bicep
L87032	Tendonitis;bicep
W34010	Test;anti-Mullerian hormone
L34004	Test;crosslaps
L41124	Test;densitometry
A33060	Test;trichomoniasis
A35010	Test;urine;metal
L57011	Therapy;laser;musculoskeletal
X99038	Torsion;ovary
A86024	Toxicity;heavy metal
T02007	Unable (to);stop eating

Should you require any assistance, please do not hesitate to contact the Medtech Helpdesk on 1800 148 165 → Option 1, or email ausupport@medtechglobal.com.