



**Medtech32 Australia
ICPC-2 PLUS Update
(October 2014)**

Release Notes



This Release Notes contain important information for all Medtech32 users. Please ensure the Release Notes are circulated amongst all your staff. We suggest these should be filed safely for future reference.

Table of Contents

For further information on this release, or any other queries regarding the ICPC-2 PLUS Update (October 2014), please contact the Medtech Helpdesk on 1800 148 165 → Option 1, or email support@medtechglobal.com.

Pre-requisites	3
Introduction	5
Installation	6
New Classification Codes	9
Inactivated Classification Codes	9

Pre-requisites

Please review the following pre-requisites and ensure they are met prior to running the Update:

- This Update can ONLY be run by CURRENT ICPC-2 PLUS Subscribers.

NOTE: For subscriptions enquiries, please contact the Family Medicine Research Centre at the University of Sydney on **02 9354 0616**.

- Ensure the person(s) who will be performing the upgrade have READ THROUGH the release notes.

IMPORTANT: This document contains valuable information that, if not read, could seriously affect the upgrade progress and/or possible down time of your network.

- Ensure you are ALWAYS logged onto Windows with ADMINISTRATIVE RIGHTS when performing ANY installation, update, or maintenance tasks.
- "User Account Control" (UAC) MUST be DISABLED as a policy across ALL Server and Client computers that are running on Windows Vista or Windows 7 or Windows 2008 or Windows 2008 R2.
- To avoid unnecessary problems during installation, upgrade, or maintenance, preferably you should ALWAYS log onto Windows in "Console Mode" – i.e. not through Remote Desktop Connection, Terminal Services, or Citrix.
- Ensure ALL Briefcasing Laptops with OUTSTANDING Briefcased data are CHECKED-IN prior to running the update.
- Ensure you have a COMPLETE backup of ALL databases located in the MT32\Data directory.

IMPORTANT: Before doing any Copy and Paste of the MT32.IB and BLOB.IB database files, you must first Stop the Interbase Server Manager. Once you have successfully copied the files to a different location, you can Restart the Interbase Server Manager.

- Ensure **ALL** activities that require access to the databases have been **STOPPED**.
 - ✓ **ALL Users** have **LOGGED OUT** of Medtech32 – including remote users and idle users in Terminal/Citrix sessions
 - ✓ **ALL Services** have been **STOPPED** – e.g. Medtech Services (for ManageMyHealth SMS and Portal), Medtech Data Transfer Service (for RxSQL Clinical Link), IIS Server (for Medtech Clinical Audit Tool), eRx Services, MediSecure Services, etc.
 - ✓ **ALL Applications** have been **STOPPED** – e.g. Medtech CDA Bridge (for NEHTA HI Service, SMD, PCEHR), eClinic SMSC Desktop, Pen Clinical Audit Tool, other Custom/Third-party Applications and Reporting Utilities not supplied by Medtech, etc.
 - ✓ **ALL Scheduled Tasks** have been **STOPPED** – e.g. Message Transfer Utility, Auto Scan from Folder, NPS RADAR Update/Submit, etc.
 - ✓ **ALL Backup/Maintenance Tasks** have been **STOPPED** – e.g. Database File Backup, Interbase Backup and Restore, Interbase Database Sweep, etc.
- Ensure the **LAST** Database you have logged in **LOCALLY** on the Medtech32 Server is the one you would like this update to run.

HINT: If unsure, simply login to the **CORRECT** Database from the Medtech32 Server, then log out and close Medtech32 before continuing.

Introduction

These release notes provide an overview on how to apply the October 2014 ICPC-2 PLUS Update for Medtech32 ICPC-2 PLUS Subscribers.

IMPORTANT NOTE

The ICPC-2 PLUS Update DOES NOT update any new ICPC-2 PLUS Groupers for data extraction from Medtech32 to Pen Clinical Audit Tool. Medtech is currently formulating the ongoing solution with Pen Computing Systems and will advise customers of the outcome.

IMPORTANT NOTE

WARNING: It is HIGHLY recommended to employ ONLY qualified system engineers when performing ANY installation and upgrade. The consequences of ruining a database during an upgrade could possibly lead to data corruptions, and as a result – data loss and systems downtime.

If in doubt, please consult with your IT technician/service provider, or contact one of the Medtech Channel Partners listed on our web site:

<http://www.medtechglobal.com/au/support-au/channel-partners-au/>

IMPORTANT NOTE FOR BRIEFCASING

Changes from this Update will NOT be applied to the Briefcasing laptops UNTIL a **COMPLETE CHECK-OUT** has been run.

IMPORTANT: ALL Briefcasing laptops with OUTSTANDING Briefcased data MUST be checked-in PRIOR TO running a complete check-out.


Alternatively, if it is not feasible to run a complete check-out soon, you MUST install this Update separately on EACH Briefcasing laptop by following the same instructions in the "Installation" section below.

Installation

The ICPC-2 PLUS Update (October 2014) must be run on the Medtech32 Server machine. The following procedures ONLY need to be run ONCE for EACH practice (or once per database if your practice has multiple databases).


NOTE: If you are uncertain which computer is the Medtech32 Server, please contact your IT technician or service provider who has performed the Medtech32 installation and/or upgrade.

1. Please visit our website at www.medtechglobal.com.
2. Select **Region: Australia** from the Region dropdown list on the top right corner of the screen.



The Australia Home Page will be displayed.

3. Select from the top menu, **Support ► Medtech32 Downloads ► Third Party Database Updates**. The Medtech32 Third Party Database Updates page will be displayed.



4. Here you will find the **October 2014 ICPC-2 PLUS Update**. Click on the title **October 2014 ICPC-2 PLUS Update** to expand the section.
5. Under the **File Attachments** section, download and **READ THROUGH** the **Release Notes**.
6. Under the **File Attachments** section, download and **RUN** the **Installer**.

7. If the File Download Security Warning dialogue box appears, select the **Run** option.
8. If the Open File Security Warning dialogue box appears, select the **Run** option.
9. The Update will then begin to run.



10. You will then be prompted to enter the **Password**. This Password has been supplied to you separately via fax.




Note: the password is not case sensitive.


NOTE: You are only eligible to receive the Password if you are currently paying for ICPC-2 PLUS subscription. If you are a current ICPC-2 PLUS subscriber but have not received the Password by fax, please contact the Medtech Helpdesk on 1800 148 165 → Option 1.

NOTE: If are not currently paying for ICPC-2 PLUS subscription but would like to re-subscribe in order to obtain this Update, please contact the Family Medicine Research Centre at the University of Sydney on 02 9354 0616.

- Enter the **Password**, and then click on the **OK** button to start the Update. The ICPC-2 PLUS Update Utility will begin to run.



- Once the Update has successfully completed, the Update Complete screen will be displayed. Click on the **OK** button to exit the Update.



New Classification Codes

The following is the list of new ICPC-2 PLUS Term Codes added by this Update:

Z62007	Notification;child abuse
R49003	Asthma action plan
D99066	Enteropathy
K76020	Aneurysm;ventricular
L81063	Haematoma;subperiosteal
R53003	Intubation;respiratory
K73015	Patent foramen ovale
D29024	Lump(s);anus
L07007	Dysfunction;temporomandibular
P17009	Problem;smoking;e-cigarette
L91023	Arthritis;toe(s)
Y45009	Advice/education;circumcision
X21018	Lesion;breast;F
Y16016	Lesion;breast;M
S99068	Disease;Grovers
H51010	Myringotomy
Y54007	Repair;penile chordee
A45041	Advice/education;test(s)
K54021	Repair;aneurysm
P24008	Impairment;communication

Inactivated Classification Codes

The following is the list of existing ICPC-2 PLUS Term Codes inactivated by this Update:

Z66003	Notification;child abuse
--------	--------------------------

Should you require any assistance, please do not hesitate to contact the Medtech Helpdesk on 1800 148 165 → Option 1, or email support@medtechglobal.com.